



## EyeMax2™

### Phototherapy Eye Masks

The soft, safe and secure solution for eye protection during neonatal phototherapy



- Adjustable band for optimum comfort and a secure fit
- Easy to use design
- Ergonomically designed eye protection pad to reduce pressure and increase light protection

The EyeMax2 has been designed to assure optimum comfort and eye protection for neonates during phototherapy treatment. The single-piece, wrap-around design is secured at two points on either side of the neonate's head and can be independently adjusted to prevent unwanted movement and allows a comfortable and secure fit.

The large eye-protection pad has been ergonomically designed to reduce pressure and increase light protection. The headband; however, allows blue light to penetrate, so that the head continues to receive phototherapy treatment.

Part Number	Size	Occipital - Frontal Circumference
1114005	Regular	32 - 38 cm
1114006	Premie	26 - 32 cm
1114007	Micro	20 - 26 cm

All EyeMax2 masks are supplied in packs of 20

Distributed by:



Part Number: 0491002  
Date: 07/12



#### Other phototherapy products available from Viamed

An alternative to the traditional phototherapy eye mask is the phototherapy light shield. These light shields are manufactured by Viamed using Amber 300™ material, and are available in a variety of standard designs and sizes. Light shields can also be designed to meet exact specifications, should the standard range not meet individual requirements.



**Viamed Limited** · 15 Station Road · Cross Hills  
Keighley · West Yorkshire · BD20 7DT · United Kingdom  
Tel: +44 (0)1535 634 542 Fax: +44 (0)1535 635 582  
Email: [info@viamed.co.uk](mailto:info@viamed.co.uk) Website: [www.viamed.co.uk](http://www.viamed.co.uk)



BS EN ISO 9001:2008  
ISO 13485:2003

FS 28344 MD 78787